



CERTIFIED UNDER:

- ANSI/ISEA Z87.1-2020 (tested in American laboratory)
- GB 32166. 1-2016

PRODUCT DESCRIPTION:

- Designed for procedures and environments where risk of exposure to blood and/or bodily fluids through splashing is possible
- Made from transparent PET (Polyethylene Terephthalate)
- Double sided anti-fog treatment
- Protection film on both sides for ultimate visibility (needs to be removed prior use)
- Equipped with foam band and headband
- Can be worn with or without glasses or goggles
- 24 cm high x 36 cm wide (9.44" x 14.17") as per required Canadian standard.

1- PACKAGING INFORMATION (100 units):

- 10 pieces per bag
- 10 bags per carton
- 100 pieces per carton
- Carton weight: 11.66 lbs (5.29 kg)
- Carton dimensions: 14.7" x 13.3" x 16.7" (37.34 x 33.78 x 42.42 cm)

2- PACKAGING INFORMATION (20 units):

- 10 pieces per bag
- 2 bags per carton
- 20 pieces per carton
- Carton weight: 5.42 lbs (2.46 kg)
- Carton dimensions: 15.38" x 14.5" x 2.5" (39.07 x 36.83 x 6.35 cm)

1 pack of 10 units	\$4.99/unit
1 carton (100 units)	\$4.49/unit
10 cartons (1000 units)	\$3.87/unit

*USD currency