

NON-CONTACT INFRARED FOREHEAD THERMOMETER

CERTIFIED UNDER:

- GB 9706.1-2007
- GB/T 14710-2009
- EN 60601-1:2006+A1:2013+A11:2011+A12:2014
- EN 60601-1-2:2015

PRODUCT DESCRIPTION:

- **Fast and accurate:** One second temperature measurement, accurate to $\pm 0.3^{\circ}\text{C}$.
- **Contactless temperature measurement:** Convenient and practical 1.18" to 1.97" (3 to 5 cm) measuring distance.
- **Multi-Purpose:** Measurable temperature of 5 to 80°C for human body and objects.
- **3 different backlight colors:** Green (normal), orange (low fever), red (high fever).

SPECIFICATIONS:

- **Size:** 3.74" x 1.65" x 6.14" (9.5 x 4.2 x 15.6 cm)
- **Weight:** 0.27 lbs (124 grams)
- **Measuring distance:** 1.18" – 1.97" (3-5 cm)
- **Temperature range:**
 - **Body:** 89.6°F- 109.2°F (32°C-42.9°C)
 - **Object:** 41°F-176°F (5°C-80°C)
- **Accuracy:**
 - 89.6°F~94.8°F: $\pm 32.5^{\circ}\text{F}$ (32°C~34.9°C: $\pm 0.3^{\circ}\text{C}$)
 - 95°F~107.6°F: $\pm 32.4^{\circ}\text{F}$ (35°C~42°C: $\pm 0.2^{\circ}\text{C}$)
 - 107.8°F~109.2°F: $\pm 32.5^{\circ}\text{F}$ (42.1°C~42.9°C: $\pm 0.3^{\circ}\text{C}$)
- **Resolution ratio:** 32.2°F (0.1°C)
- **Automatic shutdown:** 60 seconds
- **Battery:** DC3V (2x AAA Alkaline battery, number 7, not included)
- **Operating temperature:** 60.8°F-95°F (16°C-35°C)
- **Preset alarm temperature:** 100.4°F (38°C)

PACKAGING INFORMATION:

- 60 units per carton
- Carton weight: 19 lbs (8.6 kg)
- Carton dimensions: 19.29" x 14.96" x 11.42" (49 x 38 x 29 cm)



| | |
|-------------------------------------|---------------------|
| 1 unit | \$79.00/unit |
| 10 units | \$69.30/unit |
| 1 carton + (60 units/carton) | \$59.25/unit |

*USD currency