



CERTIFIED UNDER:

- ISO 9001
- EC1907-2006
- EN14683:2019+AC:2019
- YY/T 0969-2013
- EC1907-2006 (TUV)

PRODUCT DESCRIPTION:

Type II masks are primarily intended for use by healthcare professionals in an operating room or other medical setting where the same conditions are required. The two types (I & II) differ according to their bacterial filtration efficiency. Type II is further subdivided depending on whether or not the mask is splash resistant (R: splash resistant). Type II bacterial filtration efficiency (BFE) is $\geq 98\%$.

- Medical
- Non-sterilized
- With ear loops – 3 layers
- Mask dimensions: 6.89" x 3.74" (17.5 x 9.5 cm)
- Style: Flat surgical mask

PACKAGING INFORMATION:

- 50 pieces per box
- Inner box dimensions: 7.87" x 2.6" x 4.13" (20 x 6.6 x 10.5 cm)
- 40 boxes per carton
- 2000 pieces per carton
- Carton dimensions: 16.54" x 14.37" x 17.72" (42 x 36.5 x 45 cm)
- Carton weight: 21.38 lbs (9.7 kg)



| | |
|----------------------------------|--------------------|
| 1 box of 50 units | \$0.95/unit |
| 1 carton (2000 units) | \$0.93/unit |
| 10 cartons (20 000 units) | \$0.88/unit |
| 500 cartons + | \$0.78/unit |

*USD currency